





Tongue Holding Forceps, Catheter Introducing Forceps



Mouth Gags



Laryngoscope Sets, McIntosh



PAD. 23.40 – PAD. 23.55



PAD. 43.40 – PAD. 43.55



PAD. 22.00 – PAD. 23.05
PAD. 42.00 – PAD. 45.05

MCINTOSH (CONVENTIONAL)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 22.00	PAD. 23.40	Set of 3 blades No. 2,3,4
PAD. 22.03	PAD. 23.53	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 23.04	PAD. 23.54	Set of 4 blades No. 1,2,3,4
PAD. 23.05	PAD. 23.55	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps

MCINTOSH (FIBEROPTIC)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 42.00	PAD. 43.40	Set of 3 blades No. 2,3,4
PAD. 42.03	PAD. 43.50	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 43.04	PAD. 43.54	Set of 4 blades No. 1,2,3,4
PAD. 43.05	PAD. 43.55	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps

Laryngoscope Sets, Wis-Forregger, Miller



PAD. 26.20 – PAD. 26.26
PAD. 46.40 – PAD. 46.56

WIS-FOREGGER (CONVENTIONAL)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 26.00	PAD. 26.20	Set of 3 blades No. 2,3,4
PAD. 26.04	PAD. 26.24	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 26.05	PAD. 26.25	Set of 4 blades No. 1,2,3,4
PAD. 26.06	PAD. 26.26	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps

WIS-FOREGGER (FIBRE OPTIC)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 46.00	PAD. 46.40	Set of 3 blades No. 2,3,4
PAD. 46.04	PAD. 46.54	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 46.05	PAD. 46.55	Set of 4 blades No. 1,2,3,4
PAD. 46.06	PAD. 46.56	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps



PAD. 28.20 – PAD. 28.28
PAD. 48.40 – PAD. 48.56

MILLER (CONVENTIONAL)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 28.00	PAD. 28.20	Set of 3 blades No. 1,2,3
PAD. 28.04	PAD. 28.24	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 28.05	PAD. 28.25	Set of 4 blades No. 1,2,3,4
PAD. 28.06	PAD. 28.28	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps

MILLER (FIBRE OPTIC)

Wooden Box	Artificial Leather Emergency Bag	Contents
PAD. 48.00	PAD. 48.40	Set of 3 blades No. 1,2,3
PAD. 48.04	PAD. 48.54	Set of 4 blades No. 0,1,2,3 (Pediatric)
PAD. 48.05	PAD. 48.55	Set of 4 blades No. 1,2,3,4
PAD. 48.06	PAD. 48.56	Set of 5 blades, No. 1,2,3,4,5

All sets include 1 standard battery handle
and 2 spare lamps

McIntosh Blades, Conventional



PAD. 24.10
Fig. 0, 55 mm
For Premature Babies



PAD. 24.11
Fig. 1, 70 mm
For Infants



PAD. 24.12
Fig. 2, 90 mm
For Children



PAD. 24.13
Fig. 3, 110 mm
For Adults, Medium



PAD. 24.14
Fig. 4, 130 mm
For Adults, Large



PAD. 24.15
Fig. 5, 150 mm
For Adults, Extra Large

McIntosh Blades, Fibre Optic (Cold Light)



PAD. 44.10
Fig. 0, 55 mm
For Premature Babies



PAD. 44.11
Fig. 1, 70 mm
For Infants



PAD. 44.12
Fig. 2, 90 mm
For Children



PAD. 44.13
Fig. 3, 110 mm
For Adults, Medium



PAD. 44.14
Fig. 4, 130 mm
For Adults, Large

Wis-Foregger Blades, Conventional



PAD. 26.10
Fig. 0, 55 mm
For Premature Babies



PAD. 26.11
Fig. 1, 70 mm
For Infants



PAD. 26.12
Fig. 2, 90 mm
For Children



PAD. 26.13
Fig. 3, 110 mm
For Adults, Medium



PAD. 26.14
Fig. 4, 135 mm
For Adults, Large

Wis-Foregger Blades, Fibre Optic (Cold Light)



PAD. 46.10
Fig. 0, 55 mm
For Premature Babies



PAD. 46.11
Fig. 1, 70 mm
For Infants



PAD. 46.12
Fig. 2, 90 mm
For Children



PAD. 46.13
Fig. 3, 110 mm
For Adults, Medium



PAD. 46.14
Fig. 4, 135 mm
For Adults, Large

Miller Blades, Conventional



PAD. 28.30
Fig. 0, 55 mm
For Premature Babies



PAD. 28.31
Fig. 1, 80 mm
For Infants



PAD. 28.32
Fig. 2, 130 mm
For Children



PAD. 28.33
Fig. 3, 170 mm
For Adults, Medium



PAD. 28.34
Fig. 4, 180 mm
For Adults, Large

Miller Blades, Fibre Optic (Cold Light)



PAD. 48.30
Fig. 0, 55 mm
For Premature Babies



PAD. 48.31
Fig. 1, 80 mm
For Infants



PAD. 48.32
Fig. 2, 130 mm
For Children



PAD.48.33
Fig. 3, 170 mm
For Adults, Medium



PAD. 48.34
Fig. 4, 180 mm
For Adults, Large

Battery Handles & Lamps

PAD. 29.66



PAD. 29.65



FIBRE OPTIC
PAD. 30.40
ø 24 mm
11 cm / 4 1/4"



STANDARD
PAD. 31.11
ø 19 mm
13 cm / 5"



STANDARD
PAD. 32.12
ø 28mm
14 cm / 5 1/2"



FIBRE OPTIC
PAD. 32.62
ø 28mm
14 cm / 5 1/2"